



MULTI STYM TENS PROFI

MULTI STYM/TENS PROFI family of electrostimulators represents a state-of-the-art in muscle stimulation and pain relief therapy. The devices offer the user a wide range of applications in sports medicine, rehabilitation, vascular and pain relief therapy, and in wellness programs. They enable a safe and effective treatment of healthy as well as damaged muscles, increase local blood and lymph circulation, reduce the local fat, and help relieve acute and chronic pain.

Large display, a selection of programs covering a variety of indications, simple setting of treatment parameters, and small size make the MULTI STYM/TENS PROFI devices extremely user-friendly.

Although they are basically intended to be used by the professionals, their easy-to-use operation makes them a valuable tool also for home use.



CE 005

1700 224 MULTI STYM 4/220 (EMS – FES – TENS)

4 channels, 220 programs, 1000 applications intended for use in:



CE 005

- ✓ Sports medicine (increase in muscle strength and mass, increase in muscle explosive strength and muscle endurance; muscle relaxation),
- ✓ Wellness (increase in muscle tone; muscle relaxation; anti-stress therapy),
- ✓ Vascular therapy (improved local blood circulation; cramp relief; lymph drainage; relaxing massage; leg pain),
- ✓ Rehabilitation (atrophy; muscle strength; incontinence),
- ✓ Pain relief (cervical, back, lower back pain; lumboschialgia; joint pain).

Basic set: 4 x electrode cable, rechargeable battery Ni-MH 1.3 AH 7.2 V, battery charger, 4 x small self-adhesive electrode, 2 x large self-adhesive electrode, protective carrying case

Pulse shape: Symmetric rectangular biphasic; TENS; TENS modulated

Intensity: 0 – 120 mA, adjustable

Frequency: 1 – 150 Hz

Pulse duration: 30 – 450 µs

Dimensions: 145 x 97/90 x 40/35 mm

Weight: 400 g

1700 225 MULTI STYM 4/190 (EMS – FES – TENS)

4 channels, 190 programs, 750 applications intended for use in:

- ✓ Sports medicine (increase in muscle strength and mass, increase in muscle explosive strength and muscle endurance; muscle relaxation),
- ✓ Wellness (increase in muscle tone; muscle relaxation; anti-stress therapy),
- ✓ Vascular therapy (improved local blood circulation; lymph drainage; relaxing massage; leg pain),
- ✓ Rehabilitation (atrophy; muscle strength),
- ✓ Pain relief (cervical, back, lower back pain; lumboschialgia; joint pain).

Basic set: 4 x electrode cable, rechargeable battery Ni-MH 1.3 AH 7.2 V, battery charger, 4 x small self-adhesive electrode, 2 x large self-adhesive electrode, protective carrying case

Pulse shape: Symmetric rectangular biphasic; TENS; TENS modulated

Intensity: 0 – 120 mA, adjustable

Frequency: 1 – 150 Hz

Pulse duration: 30 – 450 μ s

Dimensions: 145 x 97/90 x 40/35 mm

Weight: 400 g



CE

1700 226 MULTI STYM 2/200 (EMS – FES – TENS)

2 channels, 200 programs, 800 applications intended for use in:

- ✓ Sports medicine (increase in muscle strength and mass, increase in muscle explosive strength and muscle endurance; muscle relaxation),
- ✓ Wellness (increase in muscle tone; muscle relaxation; anti-stress therapy),
- ✓ Vascular therapy (improved local blood circulation; cramp relief; lymph drainage; relaxing massage; leg pain),
- ✓ Rehabilitation (atrophy; muscle strength; incontinence),
- ✓ Pain relief (cervical, back, lower back pain; lumboschialgia; joint pain).

Basic set: 2 x electrode cable, rechargeable battery Ni-MH 1.3 AH 7.2 V, battery charger, 4 x small self-adhesive electrode, 2 x large self-adhesive electrode, protective carrying case

Pulse shape: Symmetric rectangular biphasic; TENS; TENS modulated

Intensity: 0 – 100 mA, adjustable

Frequency: 1 – 150 Hz

Pulse duration: 30 – 450 μ s

Dimensions: 145 x 97/90 x 40/35 mm

Weight: 400 g



CE

1700 227 MULTI STYM 2/110 (EMS – FES – TENS)

2 channels, 110 programs, 390 applications intended for use in:

- ✓ Sports medicine (increase in muscle strength and mass, increase in muscle explosive strength and muscle endurance; muscle relaxation),
- ✓ Wellness (increase in muscle tone; muscle relaxation),
- ✓ Vascular therapy (improved local blood circulation; lymph drainage; relaxing massage),
- ✓ Rehabilitation (atrophy; muscle strength),
- ✓ Pain relief (cervical, back, lower back pain; lumboschialgia; joint pain).

Basic set: 2 x electrode cable, rechargeable battery Ni-MH 1.3 AH 7.2 V, battery charger, 4 x small self-adhesive electrode, 2 x large self-adhesive electrode, protective carrying case

Pulse shape: Symmetric rectangular biphasic; TENS; TENS modulated

Intensity: 0 – 100 mA, adjustable

Frequency: 1 – 150 Hz

Pulse duration: 30 – 450 μ s



CE

МЕД СОЮЗ

+7 (495) 308-9009 многоканальный, www.medsoyuz.ru E-mail: ru@medsoyuz.ru

1700 228 TENS PROFI (TENS – IONOFOREZA)

2 channels, 23 programs intended for:

- ✓ Pain relief (cervical, back, lower back pain; lumboschialgia; shoulder, knee pain; tennis elbow; iontophoresis).

Basic set: 2 x electrode cable, battery 9 V, 4 x small self-adhesive electrode, protective carrying case

Pulse shape: Symmetric rectangular biphasic; TENS; TENS modulated; iontophoresis

Intensity: 0 – 60 mA, adjustable

Frequency: 1 – 150 Hz

Pulse duration: 30 – 450 μ s

Dimensions: 145 x 97/90 x 40/35 mm

Weight: 250 g



Distribution:

Iskra Medical d.o.o.,

Stegne 23, 1000 Ljubljana, Slovenia,

tel.: ++386 (0)1 513 15 00, 513 15 04,

fax: ++386 (0)1 513 15 10,

Mobil: ++386 (0) 41 622 918,

E-mail: iskra-medical@iskramedical.si,

http://www.iskramedical.si



CE